



Recombinant Human Olfactory receptor 9K2 (OR9K2)

Product Code	CSB-CF850888HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8NGE7
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MLGSKPRVHLYILPCASQQVSTMGDRGTSNHSEMTDFILAGFRVRPELHILLFL LFLFVY AMILLGNVGMMTIIMTDPRLNTPMYFFLGNLSFIDLFYSSVIEPKAMINFWSENK SISFA GCVAQLFLFALLIVTEGFLLAAMAYDRFIAICNPLLYSVQMSTRCTQLVAGSYF CGCIS SVIQTSMTFTLSFCASRAVDHFYCDSRPLQRLSCSDLFIHRMISFSLSCIILPTII VII VSYMYIVSTVLKIHSTEGHKKAFSTCSSHLGVVSVLYGAVFFMYLTPDRFPELS KVASLC YSLVTPMLNPLIYSLRNKDVQEALKKFLEKKNIL
Source	in vitro E.coli expression system
Target Names	OR9K2
Protein Names	Recommended name: Olfactory receptor 9K2 Alternative name(s): Olfactory receptor OR12-2
Expression Region	1-335
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.