



Recombinant Rat Lipid phosphate phosphatase-related protein type 2 (Lppr2)

Product Code	CSB-CF740941RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6W5G4
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	MAGGRPHLKRSFSIIPCFVFVESVLLGIVVLLAYRLEFTDTPVHTQGFFCYDSA YAKPY PGPEAASRAPPALIALVTAGPTLTILLGELARAFFPAPPSSSPVSGESTIVSGA CCRFS PPLRRLVRFLGVYSFGLFTTTIFANAGQVVTGNPTPHFLSVCRPNYTALGCPPP SPDRPG PDRFVTDQSACAGSPSLVAAARRAFPCKDAALCAYAVTYTAMYVTLVFRVKG SRLVKPSL CLALLCPAFLVGVVRVAEYRNHWSDVLGFLTGAAIATFLVTCVVHNFQSRPH SGRRLSP WEDLSQAPTMDSPLEKNPRPAGRIRHRHGSPHPSRRTVPAVAT
Source	in vitro E.coli expression system
Target Names	Plppr2
Protein Names	Recommended name: Lipid phosphate phosphatase-related protein type 2 EC=3.1.3.4 Alternative name(s): Plasticity-related gene 4 protein Short name=PRG-4
Expression Region	1-343
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.