



Recombinant Schizosaccharomyces pombe Uncharacterized mitochondrial carrier C12B10.09 (SPAC12B10.09)

Product Code	CSB-CF611969SXV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q10442
Product Type	Transmembrane Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Sequence	MIFNYVFSNSFSFLYLMVNVQCRIKEIMYWNICLYWDVFSKLIYKVCQFGRIHI FHGK NRHQLLRRTCHFTPSKRHSAMSFEEALGAGICAGLAVDLSLFPIDTLKTRLQAKG GFVKNG GFHGVYRGLGSILVGSAPGASLFFTTYENMKSRLSQSGLGLSDPQIHMCSASL GEIAACI VRVPTTEVIKQRAQASGGTLSSRNILQTILKSNNVWRDFYAGYGITAREIPFTLIQ FPIW EHLKWKRIKHSRNKNLAHEAAISGSIAGGIAAALTTPFDVVKTRIMTSQQRLSY VFTIK SIVAHEGFLALYKGIVPRVLWLSGGGAIFLGCYDVLNFMKA EGL
Source	in vitro E.coli expression system
Target Names	SPAC12B10.09
Protein Names	Recommended name: Uncharacterized mitochondrial carrier C12B10.09
Expression Region	1-345
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.