



Recombinant Danio rerio Solute carrier family 25 member 43 (slc25a43)

Product Code	CSB-CF713167DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5U3V7
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	MATVKKDARLTSSQSLMCGVGFAGIFSKTVTSPLEVVKILSQVGTFFHCKRGFLHS FVLICQ NEGLRAFWKGNMVSCLRLFPYSIAIHLATYKNIVNLHIDELGDISQWRAIVAGGL AGISAA LATYPLEVVETRLIAQNCQEPTYRGLLHLSLSVIYRNEGLQALYRGFSLTVLGAV PFSVGC YAVYINLDKLWQERHVRFTSLQNFINGCLAAGVAQTLSPFPFETVKKKMQAQL VLPHCGG VDVHFNGMADCFRQVIKNKGMALWSGLTANMVKIVPYFGLLFSCFEMCKQV CLYRNGYI ISPPSYKLPKPGVDQSLGPYELQEFKRYLRNRKSHKAQSSSIGNRW
Source	in vitro E.coli expression system
Target Names	slc25a43
Protein Names	Recommended name: Solute carrier family 25 member 43
Expression Region	1-345
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.