



Recombinant Pig P2Y purinoceptor 2 (P2RY2)

Product Code	CSB-CF733453PI
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5YA25
Product Type	Transmembrane Protein
Immunogen Species	Sus scrofa (Pig)
Sequence	MATGPDLWNGTVNGTSDGDDWGYRCRFHEDFKYVLLPLSYGVVVCVLGLSLN AGALYIFLC RLKTWNASTTYMFHLAVSDALYAASLPLLVYYYARGDHWPFFSTALCKLVRFLF YTNLYCS ILFLTICISVHRCLGVLRPLRSLRWGHARYARRVAAAVWGLVLACQAPALYFITT TAQGGR ITCHDTSAPELFSHFVAYSLVMLSVLFAAPFAVILVCYALMARRLLRPAYGTAG GLPRAK RKSVRTIAVVLAVFALCFLPFHVTRTLYYSFRTLDSLCHTLDAINMAYKITRPLAS ANSC LDPVLYFLAGQRLVRFARDAKPPTDATPTAQACRRLGLRRSHGTDTKRTEDS ASSEDSRR TEITPARGEDIRL
Source	in vitro E.coli expression system
Target Names	P2RY2
Protein Names	Recommended name: P2Y purinoceptor 2 Short name= P2Y2
Expression Region	1-373
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.