



Recombinant Inner membrane transport permease yhhJ (yhhJ)

Product Code	CSB-CF313336SZB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0AGH2
Product Type	Transmembrane Protein
Immunogen Species	Shigella flexneri
Sequence	MRHLRNIFNLGIKELRSLGDKAMLT LIVFSFTVSVYSSATVTPGSLNLAPIAIAD MDQS QLSNRIVNSFYRPWFLPPEMITADEMDAGLDAGRYTFAINIPPNFQRDVLAGR QPDIQVN VDATRMSQAFTGNGYIQNIINGEVNSFVARYRDNSEPLVSLETRMRFPNLDP AWFGGVM AIINNITMLAIVLTGSALIREREHGTVEHLLVMPITPFEIMMAKIWSMGLVVLVVS GLSL VLMVKGVLGVPIEGSIPLFMLGVALSLFATTSIGIFMGTIARSMPQLGLLVILVLL PLQM LSGGSTPRESMPQMVDIMLTMPPTTHFVSLAQAILYRGAGFEIVWPQFLTLMAI GGAFFT IALLRFRKTIGTMA
Source	in vitro E.coli expression system
Target Names	yhhJ
Protein Names	Recommended name: Inner membrane transport permease yhhJ
Expression Region	1-374
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.