



Recombinant Escherichia coli Inner membrane protein yibH (yibH)

Product Code	CSB-CF360326ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0AFV0
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli (strain K12)
Sequence	MDLLIVLTYVALAWAVFKIFRIPVNQWTLATAALGGVFLVSGLILLMNYNHPYTF TAQKA VIAIPITPQVTGIVTEVTDKNNQLIQKGEVLFKLDPVRYQARVDRLQADLMTATH NIKTL RAQLTEAQANTTQVSAERDRLFKNYQRYLKGSQAAVNPFSERDIDDARQNFL AQDALVKG SVAEQAQIQSQLDSMVNGEQSQIVSLRAQLTEAKYNLEQTVIRAPSNQYVTQV LIRPGTY AAALPLRPVMVFIPEQKRQIVAQFRQNSLLRLKPGDDAEVFNALPGQVFIGK LTSILPV VPGGSYQAQGV LQSLTVVPGTDGVLGTIELDPNDDIDALPDGIYAQVAVYSDH FSHVSVM RKVLLRMTSWMHYLYLDH
Source	in vitro E.coli expression system
Target Names	yibH
Protein Names	Recommended name: Inner membrane protein yibH
Expression Region	1-378
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.