



Recombinant Arabidopsis thaliana Uncharacterized protein At3g58460 (At3g58460)

Product Code	CSB-CF814490DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8LB17
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MRPNIVTEAGVQTRVGQWWNAIPFLTSSVVVVCGVIYLICLLTGYDTFYEVCF PSAIIIS RFQVYRFYTAIIFHGSLHVLFNMMALVPMGSELERIMGSVRLLYLTVLLATTNA VLHLL IASLAGYNPFYQYDHLMNECAIGFSGILFSMIVIETSLSGVTSRSVFGLFNVPAK LYPWI LLIVFQLLMTNVSLLGHLCGILSGFSYSYGLFNFLMPGSSFFTTIESASWMSSFI RRPKF IMCTGGNPSSYIPTYSAQNTTSSGFSTGNAWRSLSSWLPQREASNQSSDSR FPGRGRTL STARDPTAPAGETDPNLHARLLEDSSSPDRLSDATVNTVADSRQAPIANA AVL PQSQGRV AASEEQIQKLVAMGFDRTQVEVALAAADDDLTVAVEILMSQQA
Source	in vitro E.coli expression system
Target Names	RBL15
Protein Names	Recommended name: Uncharacterized protein At3g58460
Expression Region	1-403
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.