



Recombinant Bovine Zinc transporter 6 (SLC30A6)

Product Code	CSB-CF613126BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q0VC54
Product Type	Transmembrane Protein
Immunogen Species	Bos taurus (Bovine)
Sequence	MGTIHLFRKPKQRSFFGKLLQEFRLVAADRRSWKILLFGAINLTCTGFLLMWCSS TNSIAL TAYTYLTIFDLFSLITCLVSYWVMRKPSPAYSFGFERLEVLAVFASTVLAQLG ALFILK ESAERFLEQPEIHTGRLLVGTFFVALSFNLFTMLSIRNKPFAYVSEAASTSWLQE HVADLS RSLCGIIPGLSSIFLPRMNPFLIDLAFALCITYMLIEINNYFAVDTASAIALM TF GTMYPMSVYSGKVLLQTPPHVIGQLDKLIREVSTLDGVLEVRNEHFWTLGFG SLAGSVH VRIRRDANEQMVLAVHTNRLYTLVSTLTVQIFKDDWIRPALLSGPVAANVLNFT DHHVIP MPPLKGTDDSNPVTSTPTKPSSPPEFSFNTPGKNVSPVILLNTQTRPYGLGL NHGHAPY SSVLNQQLGVPGIGATQGFRTGFTNIPSRYGTTNNRIGQPRP
Source	in vitro E.coli expression system
Target Names	SLC30A6
Protein Names	Recommended name: Zinc transporter 6 Short name= ZnT-6 Alternative name(s): Solute carrier family 30 member 6
Expression Region	1-461
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.