



Recombinant Uncharacterized protein T28D9.4 (T28D9.4)

Product Code	CSB-CF611900CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q10023
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MVRWKKPKRFLCRKTQDYERIDDIENNDESDDDVILEINHGEHEDNYRNPM IRMAYIP IVIIAVIYSSFIFFFALFIEASAHTNSKCSEAHEKEKMSSGLLTPPMKHILEDRTKY KNE IKLLKTFNVHLTYSTWINSILEILALIYSFLLLLVDRQYLYPLSQVILYLFTRINFVYT D IYGRFFSPSDVVGAITTEVIYDLIRWQSESKLFKVFPTPCPLPLLLFPIFLSILQV FAN QKHTHVYETVLLIFFRIMTRFAGRRPFILSYQVLRSGIASFQSATSHDIAKTLNKV THCC LNVFASSAFFVLLMVLVDKQGLDKTPKCLAALYAFAMWMQLATVRYRHFIPILIL TVVIEL VITGLVSYQTGEFIFSHTAEVKDYVFTVMFTIIAVFRFVFIILFKVLIYKNPPHSLT TT VLLKPTNIPVTFSSIKAPQELDPTNPYYPQTVYNSNQSKNNLLPLTRD
Source	in vitro E.coli expression system
Target Names	T28D9.4
Protein Names	Recommended name: Uncharacterized protein T28D9.4
Expression Region	1-468
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.