



# Recombinant Human Transmembrane 7 superfamily member 4 (TM7SF4)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-CF023625HU  |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.   |
| <b>Uniprot No.</b>       | Q9H295  |
| <b>Product Type</b>      | Transmembrane Protein   |
| <b>Immunogen Species</b> | Homo sapiens (Human)  |
| <b>Sequence</b>          | <p>MGIWTSGETDIFLSLWEIYVSPRSPGWMDFIQHLGVCCLVALISVGLLSVAACWF<br/> LPSIIA<br/> AAASWIITCVLLCCSKHARCFILLVFLSCGLREGRNALIAAGTGIVILGHVENIFH<br/> NFKG<br/> LLDGMTCNLRAKSFSIHFPLLKKYIEAIQWIYGLATPLSVFDDLVS WNQTLAVSL<br/> FSPSH<br/> VLEAQLNDSKGEVLSVLYQMATTTEVLSSLGQKLLAFAGLSLVLLGTGLFMKRF<br/> LGPCGW<br/> KYENIYITRQFVQFDERERHQQRPCVPLNKEERRKYVIPTFWPTPKERKNLG<br/> LFFLPI<br/> LIHLCIWVLF AAVDYLLYRLIFS VSKQFQSLPGFEVHLKLHGEKQGTQDIIHDSS<br/> FNISV<br/> FEPNCIPKPKFLLSETWVPLSVILLILVMLGLLSSILMQLKILVSASFYPSVERKRI<br/> QYL HAKLLKKRSKQPLGEVKRRLSLYLTKIHFWLPVLK MIRKKQMDMASADKS</p> |
| <b>Source</b>            | in vitro E.coli expression system   |
| <b>Target Names</b>      | DCSTAMP   |
| <b>Protein Names</b>     | Recommended name: Transmembrane 7 superfamily member 4 Alternative name(s): Dendritic cell-specific transmembrane protein Short name= DC-STAMP FIND IL-4-induced protein  |
| <b>Expression Region</b> | 1-470   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |