



Recombinant Rat Metalloreductase STEAP3 (Steap3)

Product Code	CSB-CF735893RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5RKL5
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	MSGEMDKPLISRRLVDSGSLAEVPKEAPKVGILGSGDFARSLATRLVGSGFS VVGSRN PKRTAGLFPSLAQVTFQEEAVSSPEVIFVAVFREHYSSLCSLADQLAGKILVDV SNPTEK ERLQHRQSNAEYLASLFPACTVVKAFNVISAWALQAGPRDGNRQVLICGDQLE AKHTVSE MARAMGFTPLDMGSLASAREVEAIPRLLLPSWKVPTLLALGLFVCFYAYNFIRD VLQPYI RKDENKFYKMPLSVVNTTLPCVAYVLLSLVYLPGVLAAALQLRRGTKYQRFPD WLDHWLQ HRKQIGLLSFFFAMLHALYSFCLPLRRSHRYDLVNLAVKQVLANKSRLWVEEE VWRMEIY LSLGVLALGMLSLLAVTSIPSIANSNLNWKEFSFVQSTLGFVALMLSTMHTLTYG WTRAFE ENHYKFYLPPTFTLTLLLPCVILAKGLFLLPCLSHRLTKIRRGWERDGA VKFML PAGHT QGEKTSHV
Source	in vitro E.coli expression system
Target Names	Steap3
Protein Names	Recommended name: Metalloreductase STEAP3 EC= 1.16.1.- Alternative name(s): Six-transmembrane epithelial antigen of prostate 3 pHyde
Expression Region	1-488
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.