



# Recombinant *Oryza sativa* subsp. *indica* Photosystem II CP47 chlorophyll apoprotein (psbB)

<b>Product Code</b>	CSB-CF313552OFF
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P0C363
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
<b>Sequence</b>	MGLPWYRVHTVVLNDPGRLLSVHIMHTALVSGWAGSMALYELAVFDPSPVL DPMWRQGM FVIPFMTRLGITNSWGGWSISGGTVTNPGIWSYEGVAGAHIVFSGLCFLAAIWH WVYWDL EIFCDERTGKPSLDLPKIFGIHLFLAGVACFGFGAFHVTGLYGPGIWWSDPYGL TGKVQA VNPVWGAEGFDPFVPGGIASHHIAAGTLGILAGLFHLSVRPPQRLYKGLRMGN IETVLSS SIAAVFFAAFVWAGTMWYGSATTPIELFGPTRYQWDQGYFQQEIYRRVSDGLA ENLSLSE AWSKIPEKLAFYDYIGNNPAKGGLFRAGSMDNGDGIAVGWLGHPIFRDKEGR ELFVRRMP TFFETFPVVLVDEEGIVRADVPFRRAESKYSVEQVGTVVEFYGGELNGVSYSD PATVKKY ARRSQLGEIFELDRATLKSDGVFRSSPRGWFTFGHATFALLFFFGHIWHGART LFRDVFA GIDPDLDAQVEFGTFQKVGDPTRRQPV
<b>Source</b>	in vitro <i>E.coli</i> expression system
<b>Target Names</b>	psbB
<b>Protein Names</b>	Recommended name: Photosystem II CP47 chlorophyll apoprotein Alternative name(s): PSII 47 kDa protein Protein CP-47
<b>Expression Region</b>	1-508
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.