



# Recombinant Haloferax volcanii Protein translocase subunit SecD (secD)

<b>Product Code</b>	CSB-CF518657HTU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	D4GTK5
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Haloferax volcanii (strain ATCC 29605 / DSM 3757 / JCM 8879 / NBRC 14742 / NCIMB 2012 / VKM B-1768 / DS2) (Halobacterium volcanii)
<b>Sequence</b>	<p>MSTLRDNWRVIFLVAAILLSTFALFSPTMGNQSPIGEDDSATNLQFGLQLDGGT  RIRAPL  VGVTAEDVEFEGDSERVVERQVAAQIAGADAADIIVRTGAESSTVEATIENVTA  DDL SAA  LDAAGYAHGEVRDGVGTTRAETVRVLQSKINEAGLSGGTVQQVTTATGEHFI  LVEV PNR  DQSDVVDLVGERGTVQIDIYYPTGTDNGSRTYETREAVLTQADFTSIGTAQES  QTGSGAF  VPVSVRDDPAAEFQTAIQDTGLAQPGGTRCTYMEDGGRNTTEGCLLLVVNGE  VVNAFGMS  GGLADTMRAGEWAGAPSFQLQTRNTSEAQEIAINLRAGALPARLDLSGEDSG  TSSYISPS  QGSEFKFDLITGIVAVLAVAGVVFIRYGKQPVALPMIVTGLSEVYILLGFAAAIG  YPLD  LSVIAGFIAVIGTGVDLLIIIADEV MGE GSVKSRKVFQSRFRRAFVWIGAAAATTII  AMS PLAVLSLGLDQGF AIFTILGVIVGVLVTRPAYGDILRLLLTEDR</p>
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	secD
<b>Protein Names</b>	Recommended name: Protein translocase subunit SecD
<b>Expression Region</b>	1-524
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.