



Recombinant Escherichia coli ABC transporter ATP-binding protein Yojl (yojI)

Product Code	CSB-CF330301ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P33941
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli (strain K12)
Sequence	MELLVLVWRQYRWPFFISVMALSLASAALGIGLIAFINQRLIETADTSLLVLPFLG LLLL LMAVTLGSQALATTLGHHFVYRLRSEFIKRILDTHVERIEQLGSASLLAGLTSDV RNITI AFVRLPELVQGIILTIGSAAYLWMLSGKMLLVTAIWMAITWGGFVLVARVYKHM ATLRE TEDKLYTDFQTVLEGRKELTLNRERAEYVFNNLYIPDAQEYRHHIIRADTFHLSA VNWSN IMMLGAIGLVFWMANSLGWADTNVAATYSLTLLFLRTPLLSAVGALPTLLTAQV AFNKLN KFALAPFKAEFPRPQAFPNWQTLELRNVTFAYQDNAFSVGPINLTIKRGELLFLI GGNGS GKSTLAMLLTGLYQPQSGEILLDGKPVSGEQPEDYRKLFSAVFTDVWLFQQL GPEGKPA NPQLVEKWLAQLKMAHKLELSNGRIVNLKLSKGQKKRVALLLALAEERDIILLDE WAADQ DPHFRREFYQVLLPLMQEMGKTIFAISHDDHYFIHADRLLEMRNGQLSELTGE ERDAASR DAVARTA
Source	in vitro E.coli expression system
Target Names	yojI
Protein Names	Recommended name: ABC transporter ATP-binding protein Yojl
Expression Region	1-547
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.