



Recombinant Human MARVEL domain-containing protein 2 (MARVELD2)

Product Code	CSB-CF818704HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8N4S9
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	<p>MSNDGRSRNRDRRYDEVPSDLPYQDTTIRTHPTLHDSERAVSADPLPPPPLPL QPPFGPD FYSSDTEEPAIAPDLKPVRRFVPDSWKNFFRGKKKDPEWDKPVSDIRYISDGV ECSPPAS PARPNHRSPNSCKDPYGGSEGTFSSRKEADAVFPRDPYGSLDRHTQTVRTY SEKVEEYN LRYSYMKSWAGLLRILGVVELLLGAGVFACVTAYIHKDSEWYNLFGYSQPYGM GGVGGLG SMYGGYYTGPKTPFVLVAGLAWITTIILVLGMSMYRTILLDSNWWPLTEF GINVAL FILYMAAAIVVNDTNRGGLCYPLFNTPVNAVFCRVEGGQIAAMIFLVTMIVY LISAL VCLKLWRHEAARRHREYMEQQEINEPSLSSKRKMCEMATSGDRQRDSEVNF KELRTAKMK PELLSGHIPPGHIPKPIVMPDYVAKYPVIQTDDERERYKAVFQDQFSEYKELSA EVQAVL RKFDELDAVMSRLPHHSESRQEHERISRIHEEFKKKNDPTFLEKKERC DYLK NKLSHIK QRIQEYDKVMNWDVQGYG</p>
Source	in vitro E.coli expression system
Target Names	MARVELD2
Protein Names	Recommended name: MARVEL domain-containing protein 2 Alternative name(s): Tricellulin
Expression Region	1-558
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.