



Recombinant Escherichia coli Uncharacterized protein yfjW (yfjW)

Product Code	CSB-CF345572ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P52138
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli (strain K12)
Sequence	<p>MAFVSEREIVRKIFSKKIDFTILAFFYISSIFFLLCSGVLFQYFTAAFTKGNCYECS MKL DYIKQFYFSLETAWYLISAVAVFIASVFIQHRIKAYLTLAITWIVLTITDVALIH ALD N IAMNNILLNILYNLFGAILLSLFMCLSNLLFHLNKKIKHIPMILSAMIP LVSAMIIAII TAVIYLLFARQAVEIEMDISEGSDIAYVGVKDN EESFGFLNDKKTDTPTYLDVIK NGSLI YNDTQGLSGADIYIVSGCYALPDLLRNVPLDAKKSFLNVKKLEITQKLP MGMFI QGESAD VMPKAASRLSLSKQDDKFMLASSVTDSQIKFKSNNAQLMVAFAFMPIT TNGILH DYTIDI IINDKKYKIENHVAPLSRLDKNKKMKCEYQQISDLTNTYNINANYLTG FLLVLPK DDIIN YNNSPSVLLKTDFAFYKKTYQKLDKIYDDISNGKLSSLRATGISQFS INGKHL SLRPESE IIISEGSLYGLV NKSKKIKIYGTADLVFVDNKIMNLRKITYLQSKLEIFGSSIMDILK YI FGLGLLAISIKFIHSYFKNDVNENLFL</p>
Source	in vitro E.coli expression system
Target Names	yfjW
Protein Names	Recommended name: Uncharacterized protein yfjW
Expression Region	1-567
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.