



# Recombinant Mouse Cytochrome b-245 heavy chain (Cybb)

<b>Product Code</b>	CSB-CF733771MO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q61093
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Sequence</b>	<p>MGNWAVNEGLSIFVILVWLGLNVFLFINYYKVYDDGPKYNYTRKLLGSALALAR            APAACL            NFNCMLILLPVCRNLLSFLRGSSACCSTRIRRLDRNLTFHKMVAWMIALHTAI            HTIAHL            FNVEWCVNARVGISDRYSIALSDIGDNENEEYLNFAREKIKNPEGGLYVAVTRL            AGITGI            VITLCLILIITSSTKTIRRSYFEVFWYTHHLFVIFFIGLAIHGAERIVRGQTAESLEE            HN            LDICADKIEEWGKIKECPVPKFAGNPPMTWKWIVGPMFLYLGERLVRVFRSQ            QKVITKV            VTHPFKTIELQMKKKGFKMEVGQYIFVKCPKVSLEWHPFTLTSAPPEEDFFSIHI            RIVGD            WTEGLFNACGCDKQEFQDAWKLPKIAVDGPFPGTASEDVFSYEVVMLVGAGIG            VTPFASIL            KSVWYKYCDNATSLKLLKIYFYWLCRDTHAFEFADLLQLLETQMQRNNAN            FLSYNIYL            TGWDESQANHFVHHDEEKDVITGLKQKTLYGRPNWDNEFKTIASEHPNTTIG            VFLCGPE ALAETLSKQSSISNSESGPRGVHFIFNKENF</p>
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	Cybb
<b>Protein Names</b>	Recommended name: Cytochrome b-245 heavy chain EC= 1.-.- Alternative name(s): CGD91-phox Cytochrome b(558) subunit beta Short name= Cytochrome b558 subunit beta Heme-binding membrane glycoprotein gp91phox Neutrophil cyt
<b>Expression Region</b>	1-570
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.