



# Recombinant Haemophilus influenzae ATP-binding/permease protein CydD (cydD)

<b>Product Code</b>	CSB-CF338061HTA
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P45082
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Haemophilus influenzae (strain ATCC 51907 / DSM 11121 / KW20 / Rd)
<b>Sequence</b>	<p>MNKLQRQKYLQKWLRQQEPIKKLMRANIVLATLSSFILVAQTYFLATLLDKLIMQ            NVPRD            ELIPYFLGLIIGFGMRAILWAREKIGFQSGQLLRNHIRQKILDKIHLVGPATINQK            PAG            SWASIMLEQVENLHNFYARFLPQQSLSAIVPVVIFIAVFPLNWAAGLILMITAPLV            PLFM            IIVGIAAADNSQKNMDTLSRLSAQFLDRLRGLETLRRLFNRRTSEQTEHIENATEDF            RETTM            DVLKLAFLSSAVLEFFTSISIALMAVYFGFSYLGQIEFGTYNAPLTLFTGFFCLILA            PEF            YQPLRDLGTYYHDRAAGIGAADAIVDFLESDYLTVHQNEKTISLESAVEISAENL            VVLST            QGSALTKPLNFQIPANHNVALVGQSGAGKTSLMNVILGFLPYEGSLKINGQELR            ESNLAD            WRKHIAWVGQNPLLLQGTIKENLLLGDVQANDEEINQALMRSQAKEFTDKLGL            HHEIKDG            GLGISVGQAQRLAIARALLRKGDLLEPTASLDAQSENVLVQALNEASQHQ            TLMITH RIEDLKQCDQIFVMQRGEIVQQGKFTELQHEGFFAELLAQRQQDIQ</p>
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	cydD
<b>Protein Names</b>	Recommended name: ATP-binding/permease protein CydD
<b>Expression Region</b>	1-586
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.