



Recombinant Uncoordinated protein 58 (unc-58)

Product Code	CSB-CF631503CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q22271
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MFFYSPNVAPQPSSTSHRRPTLTHQHHLHVGHLDTRLSVQIPEILSRSLVHDP ERTRLVP EQLRFIIDENQELWTGSQPPTGFLELCASPVSSVFSNADARRMSQRKLVLRP TISIQED GKIIIDHVTARWEGANINSQSALLDADDGATVITDTIKDDQDDKEPKSCPQQTV KYIKIL TPHVILVSVLIGYLCLGAWILMLLETRTELLARSKKLVRLTNLMSNFTAESWKML NNAQH GVSNMDEGEWAATFREWMVRVSETVDDRRPIRRELNRPDDLNMHNKWNF PTAILYVLT LTTGGEVSVSDTDVGKVFSAFALVGIPLMFITAADIGKFLSETLLQFVSFWN RSVRKV KQWMSRIRHGRKSLQSTGGPNDTLDILGVDGTEELWFPIGAYVSCICIYCSI GSAMFI TWERTWSFIHAFHFGFNLIVTVGLGDIVVTDYIFLSLIVAFVIVGLSVVTMCVDLA STHL KAYFTRIHYFGRKRFLGMSEELKEIVALLGAMRRKKGGKVTWNDVRDFLDNE LRDRPFE PHELLMKLRFIDETSSGMSTIRHNSFQSDFFRESEYIRRVAALRPEQPAYL
Source	in vitro E.coli expression system
Target Names	unc-58
Protein Names	Recommended name: Uncoordinated protein 58
Expression Region	1-591
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.