



Recombinant *Schizosaccharomyces pombe* Uncharacterized membrane protein C977.17 (SPAC977.17)

Product Code	CSB-CF868457SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9P7T9
Product Type	Transmembrane Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Sequence	MSVPLRFSTPSSSPSASDNESVHDDGPTTELDTFNNTDVP RR VNTTKARQMR PKNTLKVA FSSPNLKGLDNTADSDSQPWLG GYLAGRLEDISGQSRRNYVDPY YEELNAGR RPNKPVWSLNGPLPHVLGNSVVEKISQKNQEARSRANSRVNSRA NSRANSSVSLAGMDGSPNWKRKMKSAVF GSRVKLNDEEAQLPRNKSSVSIAE QAASRPKVSFSLQSSRQPSIAEEQPQTQRKSSAITVEHAENAEPETPRNNVSF SRKPSIAEQDSSQDITMPPNEIIAEESLDSGSDTETLYLNYWCKIRHFFREGFAE FLGTLVLVVFVGSNLQATVTNGAGGSFESLSFAWFGCMLGVYIAGGISGG HVNPAVTISLAIFRKFPWYKVPYIYIFFQIWGAFFGGALAYGYHWSSITEFEGGKD IRTPATGGCLYTNPKPYVTWRNAFFDEFI GTAVLVGCLFAILDDTNSPPTQGMT AFIVGLLIAAIGMALGYQTSFTLN PARDLGPRMFAWWIGYGPHSFHLYHWWW TWGAWGGTIGGGIAGGLIYDLVIFTGPESPLNYPDNGFIDKKVHQITAKFEKEE EVENLEKTDSP IENN
Source	in vitro E.coli expression system
Target Names	SPAC977.17
Protein Names	Recommended name: Uncharacterized membrane protein C977.17
Expression Region	1-598
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.