



# Recombinant Arabidopsis thaliana ABC transporter G family member 23 (ABCG23)

<b>Product Code</b>	CSB-CF664265DOA
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q3E9B8
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Sequence</b>	MASCFHPSAMATSHREEDSIILFSASNSPDEFSSASSSSFSSSPLPTPNRYSLTV TNLSYT INHTPILNSVSLAAESSKILAVVGPSGTGKSTLLKIISGRVNHKALDPSSAVLMNN RKIT DYNQLRRLCGFVPQDDDLLPLLVKETLMYSAKFSLRDSTAKEREERVESLLS DLGLVLV QDSFVGEDEEDRGVSGGERKRVSIIVEMIRDPPIILLDEPTSGLDSRNSLQV VELLATM AKSKQRTVLFQPSYRILDYISDYILSRGSMVHLSLEHLEDIAKLGFIPEQ LNP IEFAMEIVESLRTFKPNSVAVVSSSMWPENNENDGIISKKEAFRVLVDVTEISYL CSRFC KIIYRTRKQLFLARTMQAVVAGLGLGSVYTRLKRDEEGVAERLGLFAFSLSFLLS STVEAL PIYLRERRVLMKESRGSYRISSYMIANTIAFVPFLFVVSLLFSIPVYWIVGLNPSI QAF SFFVLCVWLIIILMASSLVLFSAVSPDFISGNSLICTVLGAFFLFSGYFIPKEKIPK PWM FMYYVSLYRYPLESMVVNEYWSMREECFSSGNMGCLMTGEDVLKERGLDKD TRWINVGIM LAFFVFYRILCWGILLRKASKSTH
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	ABCG23
<b>Protein Names</b>	Recommended name: ABC transporter G family member 23 Short name= ABC transporter ABCG.23 Short name= AtABCG23 Alternative name(s): Probable white-brown complex homolog protein 24 Short name= AtWBC24
<b>Expression Region</b>	1-624
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.