



# Recombinant Parabacteroides distasonis Membrane protein insertase YidC (yidC)

<b>Product Code</b>	CSB-CF406396PYB
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	A6L9D2
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Parabacteroides distasonis (strain ATCC 8503 / DSM 20701 / CIP 104284 / JCM 5825 / NCTC 11152)
<b>Sequence</b>	<p>MDKNTLIGFLLIGVVLFAFSWFNRPTPEQLEAQRRYQDSIAKIEYAQQLELQKQ  ENKSAL  VEDSLENLPDSVRAQRLQQSFGVFGDAMVGTEYTTLQNDVVELRISNKGGR  CYARLKE  YDTYNDEPLVLFDEKESNFNFTLVATNRVVNTSDLYFTPVKGNPNSVIMRL  NTGEGSH  LDFTYTLKPDDYMVQYQILGTGLNGVLAPSTNALDLLWEQDIRQQEKGRKFED  RYVTLNY  KFMADDVEHLSESKSDSKQIPNRLKWIGYKDMFFSTVLISQEGFEATTLDSKAI  PEGDVL  KQFKTTASVPFDLQGKEATNLSFYFGPNKFALLKSFDKGVSAEQQLDLEKLVP  LGWGIFR  WVNQYFVIPLFDLGLKFIHNYGILILLMTIIVKILFPLTYKSYMSSAKMRVLRPQV  EEI  NAKYPGQDKAMERQKATMELYSRAGASPMMSGCLPMLLQMPILIALFMFFPSAI  ELRHQSF  LWAHDLSTYDAIFSWNKYIPIITPYFGNHISLFCLLMTITNIFYTKYNMEMTNTGQ  QQMP  GMKAMMYMMPLMFLVFFNQYASGLTYYYFISTLITIVQTLIFRYTINEDKLLAKL  EANKR KPMKKSFGMKRLEEAQRAQQETLRKQQEAKKKR</p>
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	yidC
<b>Protein Names</b>	Recommended name: Membrane protein insertase YidC Alternative name(s): Foldase YidC Membrane integrase YidC Membrane protein YidC
<b>Expression Region</b>	1-633
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.