



Recombinant Human Transmembrane protein 201 (TMEM201)

Product Code	CSB-CF023795HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5SNT2
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MEGVSALLARCPTAGLAGGLGVTACAAAGVLLYRIARRMKPTHTMVNCWFCN QDTLVPYG NRNCWDCPHCEQYNGFQENG DYNKPIPAQYLEHLNHVVSSAPSLRDPSPQPQ QWVSSQVLL CKRCNHHQTTKIKQLAAFAPREEGRYDEEVEVYRHHLEQMYKLCRPCQAAVE YYIKHQNR QLRALLLSHQFKRREADQTHAQNFSSAVKSPVQVILLRALAFLACAFLLTTALY GASGHF APGTTVPLALPPGGNGSATPDNGTTPGAEGWRQLLGLLPEHMAEKLCEAWA FGQSHQTGV VALGLLTCLLAMLLAGRIRLRRIDAFCTCLWALLLGLHLAEQHLQAASPSWLDT LKFSST SLCCLVGFTAATAVTRKATGPRRFRPRRFFPGDSAGLFPTSPSLAIPHPVGGG PASLFIP SPPSFLPLANQQLFRSPRRTSPSSLPGRLSRALSLGTIPSLTRADSGYLFSGSR PPSQVS RSGEFPVSDYFSLLSGSCPSPLPSPAPSVAGSVASSSGSLRHRRLISPARL NLKGQKL LLFPSPPGAPTTPSSSDEHSPHNGSLFTMEPPHVPRKPPLQDVKHALDLRSK LERGSAC SNRSIKKEDDSSQSSTCVVDTTTRGCSEEAATWRGRFGPSLVRGLLAVSLAA NALFTSVF LYQSLR
Source	in vitro E.coli expression system
Target Names	TMEM201
Protein Names	Recommended name: Transmembrane protein 201 Alternative name(s): Spindle-associated membrane protein 1
Expression Region	1-666
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.