



Recombinant Danio rerio Nucleoporin NDC1 (tmem48)

Product Code	CSB-CF742463DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q7SZC5
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	<p>MFSMKQNCWFIRKVVIIWRAVASIAWSVLLLLPITTAVFVLLSRFSLFHPIQWISDC TLLT ASSTIFSLMVLCAVVLITGFFNLEFYTLVPSIPCSRVALLGTVLHPLQCVHSLVYS SMGM LVMWCASVIISGRYSTLGTPCMQNESGDVLTCLNEYHLFLLLAGAFMGYSHSF LGVVKNM YYVSFQPIQQYKYPQFKGCLPMLLKCSVIQSLYSTRNFAALYFFFGYVPRAWIS STLNLP IDSSLQPLDSLTLGLLDFSLLYHLSISGTFLYFTWYLTVLIFRIYATEAYSFPVQSTF SED AERCLPKVVGEKSTLVMKFLALQDLALLSQHSPSRRQEVFSLSQPGGHPHNW NAISGECL CLLRDLTQRLVAHQDAVASNGRVKSQSASSDTRSASSSSSVLSGMEDVPETP RPTVPLRT PGSVFKSSVGGMHSSLTAPFTPDVDSPPFCSPAIRRLVGQQDPQSPWFGTVQS PHIMRRGP KLWSASTESQSNQSPASPAPAIAPSPPAANKKPSFLAQWLQNRKEQVKSFLAK RVLIVYLF NKLPEASSQALFADSQAHIWALQGLSHLVAASFSEDQFGVVQTTLPSILSSLV LLEAVD RHFKLPHASSKPARTVCSMGDSTYKTLRFALRAALKTAIYKITTTTFGEHLNAVNI STEHR KRLQQFLEFKE</p>
Source	in vitro E.coli expression system
Target Names	ndc1
Protein Names	Recommended name: Nucleoporin NDC1 Alternative name(s): Transmembrane protein 48
Expression Region	1-671
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.