



Recombinant Ceratitis capitata Protein white (W)

Product Code	CSB-CF613749DRF
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q17320
Product Type	Transmembrane Protein
Immunogen Species	Ceratitit capitata (Mediterranean fruit fly) (Tephritit capitata)
Sequence	MGQEDQEVLIRGGKATSTSAESLNNNNEQPYEQSSINQGFCKNYGTLSPSP ALTADNLT YSWYNLDVFGAVHQPGSSWKQLVNRVKGVFCNERHIPAPRKHLLKNDSGVA YPGELLAVM GSSGAGKTTLLNASAFRSSKGVQISPSTIRMLNGHPVDAKEMQARCAYVQQD DLFIGSLT AREHLIFQAMVRMPRHMTQKQKVQRVDQVIQDLSLGKCQNTLIGVPGRVKGL SGGERKRL AFASEALDPELLICDEPTSGLDSEMAHSVVQVLKKSQKGTVILTIHQPSSEL FELFD KILLMAEGRVAFLGTPGEAVDFFSYIGATCPTNYTPADFYVQVLAVVPGREVES RDRVAK ICDNFAVGKVSREMEQNFQKLVKSNFGKEDENEYTYKASWFMQFRAVLWR SWLSVLKEP LLVKVRLQTTMVAVLIGLIFLGQQLTQVGVMMNINGAIFLFLTNMTFQNSFATITV FTTE LPVFMRETRSRLYRCDTYFLGKTIAELPLFLVVPFLFTAIAIYPLIGLRPGVDHFFT ALAL VTLVANVSTSGYLISCACSSSMALSVGPPVIIPFLLFGGFFLNSGSVPVYFKW LSYLS WFRYANEGLLINQWADVKEITCTLSNTTCPSSGEVILETLNFSASDLPDFI GLALLI VGFRISAYIALTMRARRKE
Source	in vitro E.coli expression system
Target Names	W
Protein Names	Recommended name: Protein white
Expression Region	1-679
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.