



Recombinant Protein mxiA (mxiA)

Product Code	CSB-CF358070SZB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0A1I5
Product Type	Transmembrane Protein
Immunogen Species	Shigella flexneri
Sequence	<p>MIQSFLKQVSTKPELILVLMVMIIAMLIPLPTYLVDFLIGLNIVLAILVFMGSFYIE R ILSFSTFPSVLLITTLFRLALSISTSRLILVDADAGKIITTFGQFVIGDSLAVGFVIFS I VTVVQFIVITKGSERVAEVAARFSLDGMPGKQMSIDADLKAGIIDAAGAKERRSI LERES QLYGSFDGAMKFIKGDIAIAGIIIFVNLIGGISVGMSSQHGMSLSGALSTYTILTIGD GLV SQIPALLISISAGFIVTRVNGDSDNMGRNIMSQIFGNPFVLIVTSALALAIGMLPG FPPF VFFLIAVTLTALFYKVKVVEKEKSLSESDSSGYTGTFDIDNSHDSSLAMIENLDAI SSET VPLILLFAENKINANDMEGLIERIRSQFFIDYGVRLPTILYRTSNELKVDDIVLLINE VR ADSFNIYFDKVCITDENGDIDALGIPVVSTSYNERVISWVDVSYTENLTNIDAKIK SAQD EFYHQLSQALLNNINEIFGIQETKNMLDQFENRYPDLLKEVFRHVTIQRRISEVLQ RLLGE NISVRNLKLIMESLALWAPREKDVITLVEHVVRASLSRYICSKIAVSGEIKVVMLSG YIED AIRKGIRQTSGGSFLNMDIEVSDEVMETLAHALRELNAKKNFVLLVSVDIRRF VKRLID NRFKSILVISYAEIDEAYTINVLKTI</p>
Source	in vitro E.coli expression system
Target Names	mxiA
Protein Names	Recommended name: Protein mxiA Alternative name(s): Protein virH
Expression Region	1-686
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.