



Recombinant *Drosophila melanogaster* Protein white (w)

Product Code	CSB-CF320613DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P10090
Product Type	Transmembrane Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Sequence	<p>MGQEDQELLIRGGSKHPSAEHLNNGDSGAASQSCINQGFQAKNYGTLPPS PPEDSGSG SGQLAENLTYAWHNMDIFGAVNQPGSGWRQLVNRTRGLFCNERHIPAPRKHL LKNVCGVA YPGELLAVMGSSGAGKTTLLNALAFRSPQGIQVSPSGMRLNNGQPVDAKEMQ ARCAYVQQ DDLFIGSLTAREHLIFQAMVRMPRHLYRQRVARVDQVIQELSLSKCQHTIIGV PGRVKG LSGGERKRLAFASEALDPPLLICDEPTSGLDSFTAHSVVQVLKCLSQKGKTVI LTIHQV SSELFELFDKILLMAEGRVAFGLGTPSEAVDFFSYVGAQCPTNYPADFYVQVL AVVPGRE IESRDRIAKICDNFAISKVARDMEQLLATKNLEKPLEQPENGYTYKATWFMQFR AVLWRS WLSVLKEPLLKVRILIQTTMVAILIGLIFLGQQLTQVGMNINGAIFLFLTNMTFQ NVFA TINVFTSELPVFMREARSRLYRCDTYFLGKTIAELPLFLTVPLVFTAIAYPMIGLR AGVL HFFNCLALVTLVANVSTSGYLISCASSSTSMALSVGPPVIIPFLFGGFFLNSG SVPVY LKWLSYLSWFRYANEGLLINQWADVEPEISCTSSNTTCPSSGKVILETLNFA ADLPLD YVGLAILIVSFRVLAYLALRLRARRKE</p>
Source	in vitro E.coli expression system
Target Names	w
Protein Names	Recommended name: Protein white
Expression Region	1-687
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.