



Recombinant *Drosophila melanogaster* Neurogenic protein big brain (bib)

Product Code	CSB-CF326081DLU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P23645
Product Type	Transmembrane Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Sequence	MADESLHTVPLEHNIDYHIVTLFERLEAMRKDSHGGGGHGVNNRLSSTLQAPKR SMQAEIRTLEFWRSIISECLASFMYVFIVCGAAAGVGVGASVSSVLLATALASG LAMATLTQCFLHISGAHINPAVTLALCVVRSISPIRAAMYITAQCGGGIAGAALLY GVTVPGYQGNLQAAISHSAALAAWERFGVEFILTFLVVLICYFVSTDPMKKFMG NSAASIGCAYSACCFVSMPYLNPARSLGPSFVLNKWDSHWVYWFGPLVGGM ASGLVYEYIFNSRNRNLRHNKGSIDNDSSSIHSEDELNYDMDMEKPNKYQQS QGTYPRGQSNNGGGQAAGNGQHQAANMGQMPGVVANAGQGNYCQNLTY APPLSSKYDQQQEPLYGGTRSLYCRSPTLTRS NLNRSQSVYAKSNTAINRDIV PRPGPLVPAQSLYPMRTQQQQQQQQQQQQQVAPAPQSSHLQNNVQNQM QQRSESIYGMRGSMRGQQQPIQQQQQQQQQQQLQQQP NMGVQQQQMQ PPPQMMSDPQQQPQGFQPVYGTRTNPMDGNHKYDRRDPQQMYGVTGP RNRGQSAQSDDSSYGSYHGSAVTPPARHPSVEPSPPPPMLMYAPPPQPN AHPQPIRTQSERKVSAPVVVSQPAACAVTYTTSQGS AVTAQQQQQQQQQQ QQQQQQQQQMMMQQQQQHYGMLPLRPN
Source	in vitro <i>E.coli</i> expression system
Target Names	bib
Protein Names	Recommended name: Neurogenic protein big brain
Expression Region	1-696
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.