



Recombinant Dictyostelium discoideum ABC transporter B family member 5 (abcB5)

Product Code	CSB-CF855030DKK
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8T9W2
Product Type	Transmembrane Protein
Immunogen Species	Dictyostelium discoideum (Slime mold)
Sequence	<p>MYSINRIGIHGISNKILSRPKFFHSSSSSIFSLNQNLNNFKNNNNYILNGYNKITHN HHHQ TLFHRHYSNLKQEQENKQQNVGKNSETSPSSTTTTENITTATTATTKTDIKETQ NSTMSM LKTVPFKYLWPKDNNSKIRIITSVLLLLSAKVLTVPFIFKDIVDSLTTTQSEMLT LPL GLLLAYGAVKITSNGFQELRQTIFSKVAHDAIRDVSCSTFKRLHQLDLTFHLSR QTGSLS RIIDRGGRGINFLLNSILFHVVPTAFEISLVSYVMYTTLGWEYSALSLATIAAYTV FTVK VTKWRTQFRVKMNKMDNEASNKMMDSLINFETVKYFNDALEVERYHNYLKE YDKASLKT TSSLSFLNFGQALIFSLSMTAMMIMAAQGVVQGNLSVGDVLVNGLLFQISLPL NFLGTV YREIKQSLVDMDFLHLLNLPKISDNKDSKPLKLENGTIVFRDISFKYNDVSVQV LNNVS FECEGGKRIAIVGSSGSGKSTLLRLLYRFYDVSSGSIEIDGQDIRGIQLESLRKHI GVVP QDTVLFNDTIYYNIAYGPNATKEQVENAARAAHIHEVILNMKNGYDTVVGGERG LKLSSG EKQRVSIARAILKDSPVIFYDEATSSLDTEKEKLIMDALRELFKGRRTTIMIAHRLS TIVD ADEIIVLGTGGIILERNHQQLELEGKYRSMWLAQQ</p>
Source	in vitro E.coli expression system
Target Names	abcB5
Protein Names	Recommended name: ABC transporter B family member 5 Alternative name(s): ABC transporter ABCB.5
Expression Region	1-697
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.