



Recombinant Human DC-STAMP domain-containing protein 1 (DCST1)

Product Code	CSB-CF719396HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5T197
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MDIKHHQNGTRGQRRKQPHTTVQRLLTWGLPVSCSWFLWRQPGEFPVTALL LGAGAGLL AIGLFQLLVNPMNIYEEQKIMFLYSLVGLGAMGWGTSPHIRCASLLLVPKMLGK EGRLFV LGYALAAIYVGPVANLRHNLNNVIASLGCTVELQINNTRAAWRISTAPLRAMFK DLLSSK ELLRAETRNISATFEDLDAQVNSETGYTPEDTMDSGETAQGREARQAPASRL HLSTQKMY ELKTKLRCSYVVNQAILSCRRWFDRKHEQCMKHIWVPLLTHLLCLPMKFKFFC GIAKVME VWCRNRIPVEGNFGQTYDSL NQ SIRGLDGEFSANIDFKEEKQAGVLGLNTSW ERVSTEV DYVYRQEARLEWALGLLHVLLSCTFLLVLHASFSYMDSYNHDIRFDNIYISTYF CQIDDR RKKLGKRTLPLRKAEEKTVIFPCKPTIQASEMSNVVRELLETLPILLLLVVLGCL DWAL YSIFDTIRHHSFLQYSFRSSHKLEVKVGGDSMLARLLRKTIGALNTSSETVMES NNMPCL PQPVGLDARAYWRAAVPIGLLVCLCLLQAFGYRLRRVIAAFYFPKREKKRILFL YNDLLK KRAAFTKLRRRAILRRERQQKAPRHPLADILHRGCP LLRRWLCRRCVVCQAPE TPESYVC RTLDCEAVYCWSCWDDMRQRCVPCTPREELSSSAFSDSNDDTAYAG
Source	in vitro E.coli expression system
Target Names	DCST1
Protein Names	Recommended name: DC-STAMP domain-containing protein 1
Expression Region	1-706
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.