



Recombinant Probable calcium-binding mitochondrial carrier K02F3.2 (K02F3.2)

Product Code	CSB-CF633650CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q21153
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MSFDHLLTSSKRRELLQNLGIGHGSDVFSTGLFRKAECNAETAISQRSIPRANP DHLRPI FDRFATKEIKGKKLMTPEDFIRGYLGLYTEENYNKETVRLLASAADTTKDGDISF EEFCA FEALLCSPDALYLTAFELFDRNASDTISCDEFEAVIRHTQPLHDQDFDFSSEFIK RYFGA DKQRNVNYHSFCQLLHDFYEEQGIQAFKRYDKNGNGTISSLDFFQIMTTVKGH LLTDFVR HNLIAVSGGASGHKVTFPYAAFNSSLAKMELIKRVYVSTTRGNLDIEMTKEEF LHAIQS YTQVTPYEVEILFHLSELNHPGRKTLCLKDIQAIDPERLKRVSQMDRLINIKAVH HKDDR GVGTAFLESAYRFLGSGVAGACGATAVYPIDLVKTRMQNQRTSGSFVGEVMY KNSLDCFK KVVKFEGLLGLYRGLLPQIVGVAPEKAIKLTMNDYMRDKFTKDGKIPLYGEIAG GTGGM CQVVFTNPLEIVKIRLQTAGEVQQAGKKIGVFTVLKELGFLGLYKGSRACFLRDI PFSAI YFPAYAHAKLASADEDGMNSPGTLFASAFIAGVPAAGLVTPADVIKTRLQVAAR AGQTTY NGVIDCARKLIKEEGPMSLWKGTAAVCRSSPQFAVTLLTYEVLQRLFYVDFA GSRPTGS ELATTKTIQDESSTNPDHVG GYKLAAATFSGIEHKFGLFLPKFETSK
Source	in vitro E.coli expression system
Target Names	K02F3.2
Protein Names	Recommended name: Probable calcium-binding mitochondrial carrier K02F3.2
Expression Region	1-707
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.