



Recombinant Mouse Transient receptor potential cation channel subfamily V member 6 (Trpv6)

Product Code	CSB-CF852782MO
Abbreviation	Recombinant Mouse Trpv6 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q91WD2
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Mouse Transient receptor potential cation channel subfamily V member 6 (Trpv6)
Immunogen Species	Mus musculus (Mouse)
Sensitivity	Not Test
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MGPLQREDR PALGGANVAPGSSPVGVWHQPQPPKEPAFHPMGWSLPKEKG LILCLWNKFCRWFHRQESWAQSRDEQNLLQKRIWESPLLLAAKENDVQALS KLLKFEGCEVHQRGAMGETALHIAALYDNLEAMVLMEAAPELVFEPMTSELY EGQTALHIAVINQNVNLVRALLARGASVSARATGSVFHYRPHNLIYYGEHPLSF AACVGSEEIVRLIEHGADIRAQDSLGNLTVLHILILQPNKTFACQMYNLLLSYDG GDHLKSLELVPNNQGLTPFKLAGVEGNIVMFQHLMQKRKHIQWTYGPLTSTLY DLTEIDSSGDDQSLLELIVTTKKREARQILDQTPVKELVSLKWKRYGRPYFCVL GAIYVLYIICFTMCCVYRPLKPRITNRTNPRDNTLMQQKLLQEAYVTPKDDLRLV GELVSIVGAVIILLVEIPDIFRLGVTRFFGQTILGGPFHVIIITYAFMVLVTMVMRLT NVDGEVVPMSFALVLGWCNVMYFARGFQMLGPFTIMI QKMIFGDLMRFCWLM AVVILGFASAFYIIFQTEDPDELGHFYDYPMALFSTFELFLTIIDGPANYDVDLPF MYSVTYAAFAIATLLMLNLLIAMMGDTHWRVAHERDELWRAQVVATTVMRLER KLPRCLWPRSGICGREYGLGDRWFLRVEDRQDLNRQRIRRYAQAFQQQDGL YSEDLEKDSGEKLETARPF GAYLSFPTPSVSRSTSRSSSTNWERLRQGALRKD LRGIINRGLEDGEGWEYQI
Research Area	Metabolism
Source	in vitro E.coli expression system
Target Names	Trpv6
Expression Region	1-767aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.




Tag Info N-terminal 10xHis-tagged

Mol. Weight 88.9 kDa

Protein Length Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.