



Recombinant Human DC-STAMP domain-containing protein 2 (DCST2)

Product Code	CSB-CF729199HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5T1A1
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	<p>MPKVMKDVVHPLGGEEPSMARAVVRSVGGFTLGLSLATAYGLLELLVEGHSP WGCLVGTL TLAAFSLGMGFSRQVRATVLLLLPQAFSRQGRLLLLVAAFGLVLQGPCANTL RNFTRAS EAVACGAELALNQTAEVLQRAKQPLVSALNKIKAIARKTKEVADRVRKFFRSIM DGVKHI ARALRNVWQWLLHIGDVCNSELGNPYLKCARVFDDAKDSCMMVIPQAYHLCY VLMPFKLA LCGLASLVQVFCVIPKYIQPFLRQTIGTPVIQLLNRVRQEFEFNMTATHHFSVDL NASRS LSQVAMDLEAVSMKLHRVREALALMGFTTPLLLVLLYLQALFYRYCYLNWDH YDNIYIT SRFLRMEAVRSTAGLPTVLPLSAHEARRYIPPGSIFLSQWEKFFYILETFNLIRH LLLVL FLVFLDYAVFWVLDLARHQLQGEIVARSPVLVSLTVEGTGYAGNIYRDLVSAFD VLQQGN ISILSRRCLLRPSEPSTGYIVIGVMYGLCFFITLFGSYVSRLRRVICASYPSRE QERI SYLYNVLLSRRTNLLAALHRSVRRRAADQGHRS AFLVLASRCPLGPFVSHF WLHQAYCL GCGQPQDEGDMENVSCSTPGCQGLYCLTCFRLLDNTCSVCASPLSYQGDL DLELDSSDE EGPQLWLAQAQRKDPEQAWLLQQQLQEV LGRSLSMESTSESSDLDEEKGPQ QRKHGQQPL PEAHQPVSILTSPEPHRPPETSSATKGAPTPASEPSVPLSPPSLPDP SHPPP K</p>
Source	in vitro E.coli expression system
Target Names	DCST2
Protein Names	Recommended name: DC-STAMP domain-containing protein 2
Expression Region	1-773
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.