



Recombinant Human Parathyroid hormone 2 receptor (PTH2R)

Product Code	CSB-CF018990HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P49190
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	DSDGTITIEEQIVLVLKAKVQCELNITAQLQEGEGNCFPEWDGLICWPRGTVGK ISAVPC PPYIYDFNHKGVAFRHCNPNGTWDFMHSLNKWTWANYSDCLRFLQPDISIGKQ EFFERLY
Source	in vitro E.coli expression system
Target Names	PTH2R
Protein Names	Recommended name: Parathyroid hormone 2 receptor Short name= PTH2 receptor
Expression Region	27-145
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.