



Recombinant *Drosophila melanogaster* ER lumen protein retaining receptor (KdelR)

Product Code	CSB-CF527995DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O76767
Product Type	Transmembrane Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Sequence	MNIFRFAGDLSHVFAIIIIIIKIWKTRSCAGISGKSQILFAVVYLTRYLDLFTTYVS LYN SVMKVLFLATSGATVYLMYVKFKATYDHNHDSFRIEFLLVPCALLSLVINHEFTV MEVLW TFSIYLESVAILPQLFLVSRTGEAESITSHYLFALGSYRALYLLNWVYRYMVESH YDLIA IFAGVVQTVLYCDDFFLYITKVLKGGKQLPA
Source	in vitro E.coli expression system
Target Names	KdelR
Protein Names	Recommended name: ER lumen protein retaining receptor
Expression Region	1-212
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.