



Recombinant ER lumen protein retaining receptor (erd-2)

Product Code	CSB-CF343364CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P48583
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MNLFRTADVAAHAIAIVLLLLKIWKSRSCGISGRSQLLFALVFVTRYLDLFTNFF SFYN TAMKIFYLVASFGTVYLMWAKFKATYDRNNSFRIEFLVIPSMILALLINHEFIFM EVMW TFSIYLEAVAIMPQLFMLSRTGNAETITAHYLFALGSYRFLYILNWVYRYTESF FDPIS VVAGIVQTVLYADFFLYITRVIQSNRQFEMSA
Source	in vitro E.coli expression system
Target Names	erd-2
Protein Names	Recommended name: ER lumen protein retaining receptor
Expression Region	1-213
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.