



Recombinant Human ER lumen protein retaining receptor 3 (KDEL R3)

Product Code	CSB-CF012131HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O43731
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MNVFRILGDLSHLLAMILLGKIWRSKCKGISGKSQILFALVFTTRYLDLFTNFI SIYN TVMKVVFLLCAYVTVYMIYGKFRKTFDSENDTFRLEFLLVPVIGLSFLENYSFTL LEILW TFSIYLESVAILPQLFMISKGTGEAETITTHYLFFLGLYRALYLANWIRRYQTENFY DQIA VVSGVVQTIFYCDFFYLYVTKVLKGKKLSLPMP
Source	in vitro E.coli expression system
Target Names	KDEL R3
Protein Names	Recommended name: ER lumen protein retaining receptor 3 Alternative name(s): KDEL endoplasmic reticulum protein retention receptor 3 Short name= KDEL receptor 3
Expression Region	1-214
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.