



Recombinant Hydrogenase-4 component E (hyfE)

Product Code	CSB-CF365209SZB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0AEW3
Product Type	Transmembrane Protein
Immunogen Species	Shigella flexneri
Sequence	MTGSMIVNNLAGLMMLTSLFVISVKSRYRLSCGFYACQSLVLVSIFATLSCLFAAE QLLIW SASAFITKVLLVPLIMTYAARNIPQNIPEKALFGPAMMALLAALIVLLCAFVVQPV KLPM ATGLKPALAVALGHFLLGLLCIVSQRNILRQIFGYCLMENGSHLVLALLAWRAP ELVEIG IATDAIFAVIVMVLLARKIWRTHGTLVDVNNLTALKG
Source	in vitro E.coli expression system
Target Names	hyfE
Protein Names	Recommended name: Hydrogenase-4 component E EC= 1.-.-
Expression Region	1-216
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.