



Recombinant UPF0073 inner membrane protein yqfA (yqfA)

Product Code	CSB-CF354882SZB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P67156
Product Type	Transmembrane Protein
Immunogen Species	Shigella flexneri
Sequence	MVQKPLIKQGYSLAEEIANSVSHGIGLVFGIVGLVLLLQAVDLNASATAITSYSL YGG MILLFLASTLYHAIPHQRAKMWLKKFDHCAIYLLIAGTYTPFLLVGLDSPLARGL MIVIW SLALLGILFKLTIAHRFKILSLVTYLAMGWLSLVVIYEMAVKLAAGSVTLLAVGGV VYSL GVIFYVCKRIPYNHAIWHGFVLGGSVCHFLAIYLYIGQA
Source	in vitro E.coli expression system
Target Names	yqfA
Protein Names	Recommended name: UPF0073 inner membrane protein yqfA
Expression Region	1-219
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.