



Recombinant Drosophila melanogaster Sugar transporter SWEET1 (slv)

Product Code	CSB-CF752078DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q7JVE7
Product Type	Transmembrane Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Sequence	MSAVAYDSLLSTTAVISTVFQFLSGAMICRKYIQKKSTGDSSGVPFICGFLSCSF WLRYG VLTNEQSIVLVNIIGSTLFLVYTLIYYVFTVVKRACVKQFGFVLTVLVVVIVYTNRL EDQ RDRMIHVTGIVCCIVTVCFFAAPLASLLHVIRAKNSESPLPLIATSFVVSLLQWLI YGIL ISDSFIQIPNFLGCILSLLQLGLFVLYPPRSYSGHGYKLVEQAVPF
Source	in vitro E.coli expression system
Target Names	slv
Protein Names	Recommended name: Sugar transporter SWEET1 Alternative name(s): Protein saliva
Expression Region	1-226
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.