



Recombinant Human Lysoplasmalogenase (TMEM86B)

Product Code	CSB-CF822712HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8N661
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MDAGKAGQTLKTHCSAQRPDVCRWLSPFILSCCVYFCLWIPEDQLSWFAALV KCLPVLCL AGFLWVMSPSGGYTQLLQGALVCSAVGDACLIWPAAFVPGMAAFATAHLLYV WAFGFSPL QPGLLLLIILAPGPYLSVLQHLEPDMVLPVAAAYGLILMAMLWRGLAQGGSAG WGALLFT LSDGVLAWDTFQAQPLPHAHLVIMTTYAAQLLITLSALRSPVPKTD
Source	in vitro E.coli expression system
Target Names	TMEM86B
Protein Names	Recommended name: Lysoplasmalogenase EC= 3.3.2.2 EC= 3.3.2.5 Alternative name(s): Transmembrane protein 86B
Expression Region	1-226
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.