



Recombinant *Drosophila melanogaster* Gamma-secretase subunit Aph-1 (aph-1)

Product Code	CSB-CF892882DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9VQG2
Product Type	Transmembrane Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Sequence	MTLPEFFGCTFIAFGPPFALFVFTIANDPVRIIIILIAAAFFWLLSLLISLWYALIPL KE FLAFGVVFSVCFQEAFRYIIYRILRSTEQGLHAVAEDTRVTDNKHILAYVSGLGF GIISG MFALVNVLADMSGPGTMGLKGGTELFFVTSAAQALSILLHTFWSVIFFNAFDT NNYIHI GYVVFSLFVSLITLLNANELYTTLLINYLVTILTGVLAFRVAGGTSRSFRKFITC Q
Source	in vitro <i>E.coli</i> expression system
Target Names	aph-1
Protein Names	Recommended name: Gamma-secretase subunit Aph-1 Alternative name(s): Presenilin-stabilization factor
Expression Region	1-238
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.