



Recombinant UPF0073 membrane protein Mb1114c (Mb1114c)

Product Code	CSB-CF355481MVH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P67158
Product Type	Transmembrane Protein
Immunogen Species	Mycobacterium bovis (strain ATCC BAA-935 / AF2122/97)
Sequence	MSGQADTATTAEARTPAHAAHHLVEGVARVLTKPRFRGWIHVYSAGTAVLAG ASLVAVSW AVGSAKAGLTTLAYTAATITMFTVSATYHRVNWKSATARNWMKRADHSMIFVF IAGSYTP FALLALPAHDGRVLSIVWGGAIAGILLKMCWPAAPRSVGVPLYLLLGWVAVW YTATILH NAGVTALVLLFVGGALYSIGGILYAVRWPDPWPTTFGYHEFFHACTAVAAICHY IAMWVVF
Source	in vitro E.coli expression system
Target Names	BQ2027_MB1114C
Protein Names	Recommended name: UPF0073 membrane protein Mb1114c
Expression Region	1-242
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.