



Recombinant *Saccharomyces cerevisiae* Probable GDP-mannose transporter 2 (HVG1)

Product Code	CSB-CF315991SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0CE11
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Sequence	MIYTSSKSLQYLAVPIYTIFKNLTIILIA YGEVLFFGGKVTSMELTSFIMMVLSSVV ATW GDQQAIAIKASSLEDLDQELVESTIFVLNPGYLWMFTNCISSALFVLIMRKIRILT NFKD YDTMFYNNVLALPLLLVFSFIMEDWSTKNLSVNL SADS LAAMVISGLMSVGISY CSGWCV RVTSSTTYSMVGALNKLPIALAGLVFFDAPKNFLSFFSIFLGFLSGLLYAVAKQK KIQQQ KVLAATLEK
Source	in vitro <i>E.coli</i> expression system
Target Names	HVG1
Protein Names	Recommended name: Probable GDP-mannose transporter 2 Short name= GMT 2
Expression Region	1-249
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.