



Recombinant Arabidopsis thaliana Gamma-secretase subunit APH1-like (At2g31440)

Product Code	CSB-CF811869DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8L9G7
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MTVAAGIGYALVALGPSLSLFVSVISRKPFLILTVLSSTLLWLVSILLSGLWRPFL PLK ANVWWPYALLVITSVCFQEGLRFLFWKVYKRLEDVLDSFADRISRPRFLTDKL QIALAG GLGHGVAHAVFFCLSLLTPAFGPATFYVERCSKVFFLISAIILAFVTIHTFSMV IAFE GYAKGNKVDQIIVPVIHLTAGMLTLVNFASEGCVIGVPLLYLVASLTLVHCGKMW WQRLL ESRNQSSASR
Source	in vitro E.coli expression system
Target Names	At2g31440
Protein Names	Recommended name: Gamma-secretase subunit APH1-like
Expression Region	1-250
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.