



Recombinant *Drosophila melanogaster* Mannose-P-dolichol utilization defect 1 protein homolog (CG3792)

Product Code	CSB-CF894252DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9VMW8
Product Type	Transmembrane Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Sequence	MTDLIRQGALFLMSEKCYDNYFLYHNFLDVPCFKALLSKGLGLAIAGSVLVKVP QVLKI LNSKSGEGINIVGVLDLLAISFHLSYNFMHGYPFSAWGDSTFLAIQVTIAVLV LFFNG RKAQSGLFLVGYYVLMYVLNSGLTPMSVLFTIQSCNIPILLVGKLSQAYTNYQA GSTGQL SAATVIMMFAGSVARIFTSIQETGDFMIILTFIASTFANSVILGQLIYYWNKPAGV KVVD SKAKKPKTKKDD
Source	in vitro E.coli expression system
Target Names	CG3792
Protein Names	Recommended name: Mannose-P-dolichol utilization defect 1 protein homolog
Expression Region	1-252
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.