



Recombinant *Oryza sativa* subsp. *japonica*

Bidirectional sugar transporter SWEET6b

(SWEET6B)

Product Code	CSB-CF845082OFG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8W0K2
Product Type	Transmembrane Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Sequence	MISPDAARNVVGIIGNVISFGLFLSPVPTFWRICKRKDVEQFKADPYLATLLNCM LWVFY GIPIVHPNSILVVTINGIGLIVEGTYLFIFFLYSPNKKRLRLMLAVLGVELVFMLAVIL GV LLSAHTHKKRSMIVGILCVFFGSIMYFSPLTIMGKVIKTKSVEYMPFFLSLVCFLN GVCW TAYALIRFDIYVTIPNGLGAIFGAIQLILYACYRTPPKTKAAKDVEMPSVISGP GAAA TASGGSVVSVTVER
Source	in vitro E.coli expression system
Target Names	SWEET6B
Protein Names	Recommended name: Bidirectional sugar transporter SWEET6b Short name= OsSWEET6b
Expression Region	1-254
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.