



Recombinant Cytochrome c oxidase subunit 3 (COX-3)

Product Code	CSB-CF015074CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P24891
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MFHNFHILSLSSYAYNLFFASAGMLSSLVMFFKFGLYELFIFTLFSVLFISFAWG KDIAM EGLSGYHNFFVMDGFKFGVILFVFSEFMFFFCIFCTFFDAALVPVHELGETWSP FGMHLV NPFQVPLLNTIILLSSGVTVTWAHHSLLSNKSCTNSMILTCLLAAYFTGIQLMEY MEASF SIADGVFGSIFYLSTGFHGHVLCGGLFLAFNRLLLKNHFNYNHHHGLGLEFAILY WHFVD VVWLFLFVYVWWSY
Source	in vitro E.coli expression system
Target Names	ctc-3
Protein Names	Recommended name: Cytochrome c oxidase subunit 3 EC= 1.9.3.1 Alternative name(s): Cytochrome c oxidase polypeptide III
Expression Region	1-255
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.