



Recombinant Ascaris suum Cytochrome c oxidase subunit 3 (COIII)

Product Code	CSB-CF015074DOS
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P24879
Product Type	Transmembrane Protein
Immunogen Species	Ascaris suum (Pig roundworm) (<i>Ascaris lumbricoides</i>)
Sequence	MFHNFHILSLSSYPILIFCSSLGFTSSLVVFVKNGIFGGLLFCFLFSIFLVSFAWGK DIVM EGLSGYHNFFVMDGFKFGVLVFIFFSEFMFFFGIFWTFDAALVPAHDVGGVWS PIGMHLV NPFQVPLLNTIILLSSGVSVTWAHYSLLSNKGCANSLMLTCILAVYFTGIQLMEY KEASF SISDGIFGSIFYLSTGFHGVHVLFGGLFLFFNLLRLLMSHFNYNHHLGLEFAIIYW HFVD VVWLFLFVFWWSY
Source	in vitro E.coli expression system
Target Names	COIII
Protein Names	Recommended name: Cytochrome c oxidase subunit 3 EC= 1.9.3.1 Alternative name(s): Cytochrome c oxidase polypeptide III
Expression Region	1-255
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.